TrelliX[®]

Embolic Coil System

Instructions for Use

- English
- Français
- Deutsche
- Español
- Italiano
- Português
- Dansk
- Nederlands
- ελληνικά
- Polski
- Svenska
- русский
- Suomalainen

Shape Memory Medical Inc.
1321 Ridder Park Dr., Suite 10
San Jose, CA 95131, USA

Phone: +1 408 649 5175

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Phone: +31 45 303 0006

DEVICE DESCRIPTION

The SHAPE MEMORY MEDICAL Trellix Embolic Coil System consists of an implantable embolic coil attached to a delivery system called the Trellix Delivery Pusher.

- The Trellix Embolic Coil is a stretch-resistant platinum-tungsten alloy coil augmented with a self-expanding porous Shape Memory Polymer (SMP).
- The Trellix Delivery Pusher is a variable stiffness system with a radiopaque marker band located three (3) centimeters from its distal end for compatibility with microcatheters with 2 tip markers.

The Trellix Embolic Coil System is specifically designed for use with the Trellix Detachment Controller (Model: TDC-1, sold separately).

- The Trellix Detachment Controller provides the energy necessary to allow for electrolytic detachment of the Trellix Embolic Coil from the Trellix Delivery Pusher.

Upon deployment in the target lesion and exposure to an aqueous environment and body temperature, the porous Shape Memory Polymer (SMP) will slowly self-expand.

The Trellix Embolic Coil System and Trellix Detachment Controller are packaged separately as single units. They are provided sterile and are intended for single use only.

The Trellix Embolic Coil System is available in the following configurations:


- Trellix 18 Embolic Coil System, Complex Filler

The Trellix embolic coil contains a 300 Series stainless steel alloy that is generally considered safe.

CAUTION: Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

Figure 1 - Trellix Embolic Coil System

INDICATIONS FOR USE

The Trellix Embolic Coil System is intended to obstruct or occlude blood flow in vascular abnormalities of the neurovascular and peripheral vessels. Indications include:

- Intracranial aneurysms
- Other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae
- Arterial and venous embolizations in the peripheral vasculature

The device should only be used by physicians who have received appropriate training in interventional neuroradiology or interventional radiology on the use of the Trellix Embolic Coil System.

WARNINGS

- The Trellix Embolic Coil System is compatible with the Trellix Detachment Controller. DO NOT SUBSTITUTE any components or devices with components or devices from other manufacturers. Substitution could result in injury to patient and/or user.
- DO NOT use if the sterile package is open or damaged or if the system is damaged.
- DO NOT use after the "Use-by" date specified on the package.
- DO NOT attempt to re-sterilize or otherwise re-use.
- DO NOT advance the delivery pusher with excessive force.
- To reduce the risk of coil migration, the diameter of the first and second coil should never be less than the width of the ostium.
- Advancing the delivery pusher beyond the microcatheter tip once the coil has been detached involves risk of aneurysm or vessel perforation.
- DO NOT rotate the delivery pusher during or after delivery of the coil. Rotating the delivery pusher may result in a stretched or premature detachment of the coil from the delivery pusher, which could result in coil migration.
- Damaged delivery pushers may cause vessel injury or unpredictable distal tip response during coil deployment.
- The fluoro-saver marker is designed for use with a rotating hemostatic valve (RHV). If used without an RHV, the distal end of the coil may be beyond the alignment marker when the fluoro-saver marker reaches the microcatheter hub.
- DO NOT overtighten the RHV around the introducer sheath. Excessive tightening may cause damage to the device.
- The fluoro-saver marker is designed for use with a microcatheter length of 150cm or greater. If used with a microcatheter shorter than 150cm, the distal end of the coil may be beyond the alignment marker when the fluoro-saver marker reaches the microcatheter hub.
- Microcatheter tip movement could cause the aneurysm or vessel to rupture.
- DO NOT place the Trellix delivery pusher on a bare metallic surface.
- If implant fails to detach, to minimize disturbance to the coil placement and delivery pusher, only press the "Detach Button" again to reinitiate. DO NOT remove or modify delivery pusher position.

CONTRAINDICATIONS

None known.

PRECAUTIONS

- High quality, digital subtraction fluoroscopic road mapping is MANDATORY to achieve safe placement of the Trellix Embolic Coil System.
 - Take all necessary precautions to limit fluoroscopic radiation doses to patients and clinical operators by using sufficient shielding, reducing fluoroscopy times, and modifying x-ray technical factors whenever possible.
 - It is recommended to deploy the Trellix Embolic Coil System to the target site within ten (10) minutes of entering the microcatheter (*working time*). A *working time* greater than ten (10) minutes (i.e. exposure within an aqueous environment) may result in increased delivery friction or the device may not be able to retract into the microcatheter without risk to damage to the device.
 - The Trellix 18 embolic coil includes an outer layer of shape memory polymer (SMP). The primary coil diameter will slowly increase to a maximum 0.030" (0.76mm) upon exposure to an aqueous environment and body temperature. This may increase the delivery friction to deliver the implant to the target site.
 - In order to achieve optimal performance of the Trellix Embolic Coil System and to reduce the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate flush solution be maintained between:
 - The femoral sheath and the guiding catheter
 - The microcatheter and the guiding catheter
 - The microcatheter and the Trellix delivery pusher
- Continuous flush also reduces the potential for thrombus formation on, and crystallization of infusate around, the detachment zone of the Trellix delivery pusher.

- Advance and retract the Trellix embolic coil slowly and smoothly, especially in tortuous anatomy. Remove the entire device if excessive friction is noted. If friction is noted in a second device, carefully examine both coil and the microcatheter for possible damage such a microcatheter shaft buckling or kinking.
- Due to the delicate nature of the Trellix embolic coils, the tortuous vascular pathways that lead to certain aneurysms and vessels, and the varying morphologies of intracranial aneurysms, a coil may occasionally stretch while being maneuvered. Stretching is a precursor to potential coil breakage and migration.
- If repositioning of the Trellix embolic coil is necessary, take special care to retract the coil under fluoroscopy in a one-to-one motion with the delivery pusher. If the coil does not move in a one-to-one motion with the delivery pusher, or if repositioning is difficult, the coil may have become stretched and could possibly migrate or break. Gently remove both the coil and microcatheter and replace with new devices.
- Avoid placing the microcatheter within the aneurysm sac where the aneurysm wall may support the microcatheter positioning. This may increase the resistance required to deploy the implant into the aneurysm as the implant exits the catheter. Repositioning the distal tip of the microcatheter or placing the distal tip at the ostium of the aneurysm may be required.
- If resistance is encountered while withdrawing a coil that is at an acute angle relative to the microcatheter tip, it is possible to avoid coil stretching or breaking by carefully repositioning the distal tip of the catheter at, or slightly inside, the ostium of the aneurysm. By doing so, the aneurysm and artery act to funnel the coil back into the microcatheter.
- Multiple embolization procedures may be required to achieve the desired occlusion of some aneurysms or vessels.
- Probability of occurrence of fluoroscopic complications may increase as procedure time and number of procedures increase.
- Increased detachment times may occur when:
 - Other embolic agents are present
 - Delivery pusher and microcatheter markers are not properly aligned
 - Thrombus is present on the coil detachment zone

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur during or after a procedure include, but are not limited to:

- Allergic reaction
- Aneurysm perforation/rupture
- Arrhythmia
- Coil migration or misplacement
- Death
- Dissection
- Edema
- Embolus
- Headache
- Hemorrhage
- Hypertension
- Hypotension
- Incomplete aneurysm filling
- Infection
- Ischemia
- Neurological/intracranial sequelae
- Parent artery occlusion
- Post-embolization syndrome (fever, increased blood cell count, discomfort)
- Premature or difficult coil detachment
- Revascularization
- Stroke/TIA
- Trauma or damage
- Vasospasm
- Vessel perforation/rupture
- Vessel thrombosis
- Other procedural complications including but not limited to:
 - Access site complications
 - Access site hematoma
 - Anesthetic and contrast media risks
 - Risks associated with excessive fluoroscopic radiation such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia.

PACKAGING AND STORAGE

The Trellix Embolic Coil System is packaged in a sterile, desiccated pouch and inside a protective plastic dispenser hoop.

The Trellix Detachment Controller is packaged in a sterile pouch.

- Inspect the package and system prior to use. DO NOT use if the sterile package is open or damaged or if the system is damaged.
- Device is intended for SINGLE PATIENT USE ONLY. DO NOT attempt to re-process, re-sterilize, clean, or re-use the device. Improper sterilization and re-use can cause malfunction to the device and injury to the patient.
- DO NOT use after the "Use-by" date specified on the package.
- DO NOT use if labeling is incomplete or illegible.

Store in a cool, dark and dry place.

MR CONDITIONAL

Non-clinical testing and MRI simulations were performed to evaluate the entire family (i.e., available diameters and lengths) of the Trellix Embolic Coil System. Non-clinical testing demonstrated that the entire family of these embolization coils is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the Trellix Embolic Coil System is expected to produce a maximum temperature rise of 2.0°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the Trellix Embolic Coil System extends approximately 3-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. In addition, the use of a magnetic resonance angiography (MRA) pulse sequence to scan an aneurysm model (i.e., 22% estimated packing density of embolization coils) at 3-T, using an echo time of 3.5-ms, generated a signal void with a linear extension beyond the embolization coil mass of 6-mm.

MATERIALS RECOMMENDED

- Trellix Detachment Controller (Model: TDC-1)
- Femoral Sheath
- Guiding Catheter (5F-7F)
- Microcatheter with 2 tip markers (located 3cm apart), *minimum length: 150cm*
- Guidewire compatible with microcatheter
- Three (3) Continuous Saline (or heparinized saline) Flush set-up with pressure bags
 - One (1) flush system for the femoral sheath
 - One (1) flush system for the guiding catheter
 - One (1) flush system for the microcatheter
- Two (2) Rotating Hemostatic Valves (RHV)
- One (1) 3-Way Stopcock
- One (1) 1-Way Valve
- One (1) uncoated sterile 20- or 22-gauge stainless steel hypodermic needle (to provide ground for the detachment controller)

Note: A backup Trellix Detachment Controller is required for all procedures (each controller will allow for a minimum of twenty (20) detachment cycles).

DEVICE SELECTION

MICROCATHETER SELECTION

Proper selection of the appropriately sized microcatheter is required to avoid damage to the device and to minimize potential complications.

Microcatheter selection is also determined by the physician and is predicated by the location of the aneurysm, patient safety, and physician preference. To ensure proper placement and detachment of the device, it is recommended the selected microcatheter contain two (2) radiopaque tip markers located three (3) cm apart and is a minimum 150cm in length.

The SHAPE MEMORY MEDICAL Trellix 18 Embolic Coil System is compatible with microcatheters with inner lumen diameter of 0.021" (0.5mm).

CAUTION: The Trellix 18 embolic coil includes an outer layer of crimped shape memory polymer (SMP) approximately .015" (0.38mm) in diameter. The primary coil diameter will slowly increase to a maximum 0.030" (0.76mm) upon exposure to an aqueous environment and body temperature. This may increase the delivery friction to deliver the implant to the target site.

COIL SIZE SELECTION

Device selection is at the discretion of the physician. In order to choose the optimum device for any given lesion, examine the pre-embolization angiograms. Appropriate device size should be chosen based upon pre-embolization angiographic assessment of the diameter of the vessel, aneurysm and/or ostium. When accessing aneurysms, the diameter of the first or second coils should never be less than the width of the ostium, given the potential for these coils to migrate.

DIRECTIONS FOR USE - PREPARATION

CONTINUOUS SALINE FLUSH SET-UP

To achieve optimal performance of the Trellix Embolic Coil System, it is important that a continuous infusion of an appropriate flush solution be maintained. Figure 2 illustrates the connections necessary for the Trellix Embolic Coil System including a typical continuous saline flush set-up with pressure bag for catheter systems.

1. Attach a rotating hemostatic valve (RHV) to the hub of the guiding catheter.
 - Connect a 3-way stopcock to the RHV sidearm and then connect a line for continuous infusion of appropriate solution
2. Attach a second RHV to the hub of the microcatheter.
 - Connect a 1-way stopcock to the RHV sidearm and then connect a line for continuous flushing of appropriate solution
 - Open stopcock to allow flush through microcatheter with sterile flush solution
3. Check to ensure that all fittings are secure and that no air is being introduced into the system during active flushing.

LESION ACCESS

4. Carefully catheterize the lesion to be treated. The access system should include a guiding catheter of sufficient inner diameter (ID) to accept a microcatheter and to permit adequate contrast infusion around the microcatheter for fluoroscopic road mapping.
5. Measure the size of the lesion to be treated and select an appropriately sized Trellix Embolic Coil System.

DIRECTIONS FOR USE - TRELIX PREPARATION & PLACEMENT

TRELIX DETACHMENT CONTROLLER PREPARATION

1. Remove the Trellix Detachment Controller from its protective packaging. Place the detachment controller in the sterile field.
CAUTION: DO NOT use any power source other than the SHAPE MEMORY MEDICAL Trellix Detachment Controller to detach the SHAPE MEMORY MEDICAL Trellix Embolic Coil System. The Trellix Detachment Controller is intended to be single patient use only. DO NOT attempt to re-sterilize or otherwise re-use.
2. Insert sterile, uncoated hypodermic needle size 20 or 22 gauge (0.9mm or 0.7mm) at the patient's groin.
3. Clip the return electrode (grounding) cable (Blue cable) on the Trellix Detachment Controller onto the hypodermic needle.

TRELIX EMBOLIC COIL SYSTEM PREPARATION

4. Slowly and simultaneously remove the Trellix embolic coil with delivery pusher and introducer sheath from the dispenser hoop. Inspect proximal section of delivery pusher for irregularities. If irregularities exist, replace with a new Trellix Embolic Coil System.
5. Hold the device just distal to the shrink-lock and pull the shrink-lock proximally to expose the tab on the introducer sheath (See Figure 3).

Figure 3a - Shrink-Lock Removal


Figure 3b - Shrink-Lock Removed

6. Slowly advance the Trellix embolic coil out of the introducer sheath and inspect the implant for any irregularities or damage. If irregularities exist, replace with a new Trellix Embolic Coil System.
7. While holding the introducer vertically, gently retract the coil back into the introducer sheath about 1-2 cm.


Figure 2 - Recommended Setup Diagram

INTRODUCTION & TRANSFER OF IMPLANT

- Open the RHV on the microcatheter just enough to accept the introducer sheath.
- Insert the distal end of the introducer sheath through the RHV and into the hub of the microcatheter until the sheath is firmly seated and a minimal gap is present (See Figure 4).

Figure 4 - Introducer Alignment with Microcatheter Hub


- Close the RHV lightly around the introducer sheath to secure the RHV to the microcatheter and to prevent any back flow.

CAUTION: DO NOT overtighten the RHV around the introducer sheath. Excessive tightening may cause damage to the device.
- With the shrink-lock removed, begin transferring the Trellix embolic coil to the microcatheter by advancing the delivery pusher in a smooth, continuous manner. Use caution to avoid catching the coil on the junction between the introducer sheath and the hub of the microcatheter.

CAUTION: The shape memory polymer (SMP) has a *working time* of ten (10) minutes upon exposure to an aqueous environment and body temperature. Manipulating the device longer than ten (10) minutes within the microcatheter may risk damage to the Trellix embolic coil.
- Continue to advance until the proximal end of the introducer sheath meets the proximal end of the Trellix delivery pusher.
 - Loosen the RHV and retract the introducer sheath just out of the RHV
 - Close the RHV around the delivery pusher
 - Slide the introducer sheath completely off of the delivery pusher. Use care not to damage the delivery pusher
- Carefully advance the Trellix Embolic Coil System until the fluoro-saver markers on the proximal end of the delivery pusher approaches the RHV on the hub of the microcatheter (See Figure 5). At this time, fluoroscopic guidance must be initiated.


Figure 5 - Fluoro-saver Marker Alignment

- Under fluoroscopic guidance, slowly advance the Trellix embolic coil out of the tip of the microcatheter. Continue to advance the coil into the lesion until optimal deployment is achieved. Reposition if necessary. If the coil size is inappropriate, remove and replace with appropriately sized coil.

CAUTION: DO NOT rotate the delivery pusher during or after delivery of the coil to the lesion site. Rotating the delivery pusher may result in stretched coil or premature detachment of the coil from the delivery pusher which could result in coil migration.
- Continue to advance the coil until the radiopaque (RO) proximal marker on the delivery pusher is exactly distal to the proximal marker on the microcatheter with 2 tip markers (See Figure 6).

Figure 6 - Radiopaque Marker Alignment with Microcatheter


- Tighten the RHV to prevent movement of the delivery pusher.
- Verify that the distal shaft of the microcatheter is not under stress before coil detachment. Axial compression or tension forces could be stored in the microcatheter causing the tip to move during coil delivery.

CAUTION: Microcatheter tip movement could cause the aneurysm or vessel to rupture.

DETACHMENT OF COIL

- Verify the hypodermic needle is in place at the patient's groin and the blue cable connector of the Trellix Detachment Controller is connected to the hypodermic (grounding) needle.
- Verify that the RHV is firmly locked around the Trellix delivery pusher prior to insertion.
- Prior to insertion, wipe clean and dry the exposed proximal section of the delivery pusher.
- To insert, slide the Trellix Detachment Controller over the proximal connector of the Trellix delivery pusher firmly into the funnel section of the Trellix Detachment Controller (See Figure 7).

Figure 7 - Detachment Controller Insertion to Delivery Pusher

- When the Trellix Detachment Controller is properly connected to the Trellix Coil System, a single short audible beep will sound and a solid green light ("System Ready") will be present on the Trellix Detachment Controller indicating the system is ready for detachment. If "System Ready" light is flashing green, check to make sure that the connection has been made and return electrode (grounding) cable is connected correctly.

If "System Ready" light is blinking and no audible beep was present, verify return electrode (grounding) cable is properly connected.

If "Low Battery" indicator is illuminated amber, replace the Trellix Detachment Controller.


For all indicator descriptions and any troubleshooting, refer to the Instructions for Use included with the Trellix Detachment Controller (Model: TDC-1).
- Prior to initiating detachment:
 - Verify visually that the continuous flush solution is infusing normally
 - Verify under fluoroscopy that the proximal marker of the delivery pusher is exactly distal to the proximal marker of the microcatheter with 2 tip markers (See Figure 6).
 - Verify the coil position is stable at the lesion site
- Press and hold momentarily the "Detach Button" on the Trellix Detachment Controller to initiate detachment. A single short audible beep will sound and the "System Ready" light and "Progress Indicators" will flash green indicating detachment. Once initiated, release the button.
- During the detachment cycle, the "Progress Indicators" will progressively light up (up to a maximum duration of 120 seconds).
- At the end of a successful detachment cycle:
 - Current flow will stop
 - Three (3) short audible beeps will sound
 - All four (4) Progress Indicator lights will be a solid green light
 - System Ready light will be off
- If detachment is NOT successful due to a time-out (greater than 120 seconds):
 - Current flow will stop
 - One (1) long audible beep will sound
 - All four (4) Progress Indicator lights will flash amber
 - System Ready light will remain solid green

To reinitiate, repeat and press the "Detach Button" to initiate another cycle of detachment.

CAUTION: To minimize disturbance to the coil placement and delivery pusher, only press the "Detach Button" again to reinitiate. DO NOT remove or modify delivery pusher position.
- Once the detachment of the coil has been signaled, VERIFY UNDER FLUOROSCOPY THAT THE COIL HAS DETACHED. Slowly pull back on the delivery pusher while watching fluoroscopy to make sure the coil does not move.
- In the unlikely event the coil moves, perform the following:
 - Advance the delivery pusher to re-establish the RO marker and microcatheter marker alignment
 - Verify that the RHV is firmly locked around the Trellix delivery pusher before removing and inserting into the Trellix Detachment Controller to ensure that the coil does not move during the connection process
 - Resume current flow by pressing and holding momentarily the "Detach Button" on the Trellix Detachment Controller
 - Verify coil detachment under fluoroscopy as above. Repeat as necessary.
- Once coil detachment has been detected and fluoroscopically confirmed, remove the delivery pusher from the Trellix Detachment Controller and slowly withdraw the wire from the microcatheter.

CAUTION: DO NOT advance the delivery pusher once the coil has been detached as it may increase risk of aneurysm or vessel perforation.
- Repeat steps if additional coil placement is required.

TRELLIX DETACHMENT CONTROLLER
BASIC CONTROLLER INDICATOR DIAGRAM

Light Indicators

- System Ready (Solid Green) - Ready to initiate detachment
- System Ready (Blinking Green) - Check return electrode (grounding) cable connections
- Low Battery (Amber) - Replace controller

Audio Indicators

- Single Short Beep - Confirms correct action
- Three Short Beeps - Successful detachment
- Single Long Beep - Error alert

For full specifications, operation instructions, and troubleshooting, refer to the instructions for use (IFU) included with the Trellix Detachment Controller.

PRODUCT DISPOSAL INFORMATION

After the device usage, the Trellix Embolic Coil System should be disposed of according to normal hospital practices (e.g., biohazardous materials should be properly handled and disposed of in appropriate containers). Infected medical devices are not regulated under the WEEE Directive but rather must be disposed of as medical waste.

WARRANTY

Shape Memory Medical Inc. warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Shape Memory Medical Inc.'s obligation under this warranty is limited to replacing or repairing at its option, this product if returned within the warranty period to Shape Memory Medical Inc. and after confirmed to be defective by the manufacturer.

Except as expressly provided in this warranty, Shape Memory Medical Inc. disclaims any representation or warranty of any kind, express or implied, including any warranty as to merchantability or fitness for a particular purpose.

See the Terms and Conditions of Sale for further information.

SYMBOLS

The following symbols are used within product labeling:

	European Conformity		Do Not Resterilize
	Refer to Instructions for Use		Caution
	Manufacturer		MR Conditional
	Use-by Date		Keep Dry
	Do not use if package is damaged		Storage, Temperature
	Lot Number		Storage, Humidity
	Catalog Number		Non-Pyrogenic
	Sterilized by Radiation (Electron Beam)		Authorized European Representative
	Single Use Only. Do not Re-use		
	Federal law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).		

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