 **shape** memory medical[®]

TrelliX[®]

Detachment Controller (Model: TDC-1)

Instructions for Use

- English
- Français
- Deutsche
- Español
- Italiano
- Português
- Dansk
- Nederlands
- ελληνικά
- Polski
- Svenska
- русский
- Suomalainen

Shape Memory Medical Inc.
1321 Ridder Park Dr., Suite 10
San Jose, CA 95131, USA

Phone: +1 408 649 5175

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Phone: +31 45 303 0006

DEVICE DESCRIPTION

The SHAPE MEMORY MEDICAL Trellix Detachment Controller consists of a sterile, handheld, single patient use device designed for use with the SHAPE MEMORY MEDICAL Trellix Embolic Coil System.

The Trellix Detachment Controller provides the energy necessary to allow for electrolytic detachment of the Trellix embolic coil from the Trellix delivery pusher.

The Trellix Embolic Coil System and Trellix Detachment Controller are packaged separately as single units. They are provided sterile and are intended for single patient use only. The Trellix Detachment controller is preloaded with non-replaceable AAAA batteries and is intended for use in one clinical case, for a minimum of 20 coil detachments.

The Trellix Detachment Controller consists of a handheld detachment unit with a detachment button, a funnel for insertion of the Trellix delivery pusher, two (2) Status LEDs, four (4) Progress LEDs, and a return electrode (grounding) cable terminated with a clip. The Trellix Detachment Controller must be used in combination with a 20- or 22-gauge uncoated stainless steel hypodermic needle (not supplied) for grounding.

CAUTION: Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

INDICATIONS FOR USE

The Trellix Detachment Controller is intended for use with the Trellix Embolic Coil System which is intended to obstruct or occlude blood flow in vascular abnormalities of the neurovascular and peripheral vessels. Indications include:

- Intracranial aneurysms
- Other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae
- Arterial and venous embolizations in the peripheral vasculature

PRECAUTIONS

- Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Prior to beginning a procedure, confirm that there are enough Trellix Detachment Controllers available to complete the anticipated number of coil detachments. In addition, one (1) extra Trellix Detachment Controller is required as backup for all procedures.
- Verify that the Trellix Detachment Controller to be used is within its labeled shelf life.
- Batteries are preloaded into the Trellix Detachment Controller. Do not attempt to disassemble the unit and replace the batteries.
- Increased detachment times may occur when:
 - Other embolic agents are present
 - Detachment zone is not properly positioned for detachment
 - Trellix delivery pusher and microcatheter markers are not properly aligned
 - Thrombus is present on the detachment zone
 - The Trellix Detachment Controller's return electrode (grounding) cable is not attached firmly to the hypodermic needle, or the hypodermic needle is not placed properly within the tissue.
- Never insert a Trellix delivery pusher into the Trellix Detachment Controller's funnel unless it is free from fluids and debris, as this may permanently render the Trellix Detachment Controller unusable for future detachments or impact the system to make proper electrical contact during detachment.
- Hypodermic (grounding) needle must be uncoated.
- If the Trellix Detachment Controller is set on the table during a detachment cycle, take care to gently place it in a stable position so that it does not slide off of the delivery pusher during detachment. If the Trellix Detachment Controller is held, take care to gently hold it in a stable position so that it does not slide off of the delivery pusher during detachment.
- Because coils are not always detached following completion of a cycle, ALWAYS verify for detachment under fluoroscopy prior to removing the Trellix delivery pusher. Always confirm that the coil has detached by gently and slowly pulling back on the Trellix delivery pusher while monitoring the fluoro image that there is no coil movement. In the unlikely event the coil moves (indicating attachment to the delivery pusher), check and adjust the flush system, flush the system to clear any contrast that may exist around the detachment zone, verify the system is properly grounded, realign the delivery pusher with the microcatheter, tighten the RHV, and repeat the detachment procedure.
- If detachment times are increasing and/or the patient is experiencing pain at the needle insertion site, replace the needle with a new needle at a new insertion site.
- Portable and mobile RF communications equipment may affect the performance of the Trellix Detachment Controller.
- If implant fails to detach, to minimize disturbance to the coil placement and delivery pusher, only press the "Detach Button" again to reinitiate. DO NOT remove or modify delivery pusher position.

CONTRAINDICATIONS

None known.

POTENTIAL ADVERSE EVENTS

Refer to the Instructions for Use for the Trellix Embolic Coil System for potential adverse events associated with the coiling procedure.

Although no adverse events are expected with the Trellix Detachment Controller, failure to detach a coil may result in a longer procedure time, or additional interventional procedures.

Please notify Shape Memory Medical Inc. immediately if the Trellix Detachment Controller malfunctions or if patient injury is suspected as a result of using this system.

WARNINGS

This device should only be used by physicians trained in interventional neuroradiology or interventional radiology on the use of this device as established by Shape Memory Medical Inc.

- The Trellix Detachment Controller can only be used with the Trellix Embolic Coil System. Use of the Trellix Detachment Controller with another manufacturer's device may lead to injury to the patient or user.
- DO NOT attempt to use the Trellix Detachment Controller with another manufacturer's devices due to potential incompatibility. Safety and performance have NOT been established with another manufacturer's devices and is not recommended.
- Inspect the package and system prior to use. DO NOT use if the sterile package is open or damaged or if the system is damaged.
- The Trellix Detachment Controller is intended for SINGLE PATIENT USE ONLY. DO NOT attempt to re-process, re-sterilize, clean, or re-use the device. Improper sterilization and re-use can cause malfunction of the device and injury to the patient.
- DO NOT use after the "Use-by" date specified on the package.
- DO NOT use if labeling is incomplete or illegible.
- The Trellix Detachment Controller should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, the Trellix Detachment Controller should be observed to verify normal operation in the configuration in which it will be used.
- There are no cables, transducers or other accessories that are likely to affect compliance of this equipment with the appropriate electro-magnetic standards.
- When handling the controller during device operation, always use double-gloved hands in order to provide necessary insulation from any malfunctions that may cause heat generation.
- Advancing the Trellix Detachment Controller onto the Trellix delivery pusher with excessive force may cause a kink in the delivery pusher or displacement of the coil from the target site, which could lead to failed detachment, suboptimal coil position, vessel/aneurysm perforation, or damage to the Trellix Detachment Controller.
- DO NOT advance the delivery pusher once the coil has been detached as it may increase risk of aneurysm or vessel perforation.
- After use, the Trellix Detachment Controller and its packaging should be handled as biohazardous waste and disposed of properly, in compliance with hospital, local, and governmental laws and policies.
- No modification of the device is allowed.

PACKAGING AND STORAGE

The Trellix Detachment Controller is packaged in a sterile pouch and inside a protective shelf carton. The Trellix Detachment Controller is sterilized with Ethylene Oxide (EO).

No maintenance or calibration is required.

CONFORMANCE TO STANDARDS

Safety Information: The Trellix Detachment Controller complies with the requirements of:

- ANSI/AAMI ES60601-1:(2005)+A1(2012)
- CAN/CSA C22.2 No. 60601-1:2014 Ed.3
- IEC 60601-1:2005 Ed.3+A1(2012)
- EN 60601-1 Ed. 3.1 (2013)
- IEC 60601-1-6 (2013)
- IEC 62366-1 (2015)
- IEC 60601-1-2 Ed. 4.0 (2014)

If any interference is suspected while using the Trellix Detachment Controller, relocating or repositioning the device may reduce or eliminate the interference.

MATERIALS RECOMMENDED


- One (1) uncoated sterile 20- or 22- gauge stainless steel hypodermic needle (to provide ground for the detachment controller)

Note:

A backup Trellix Detachment Controller is required for all procedures (each controller will allow for a minimum of twenty (20) detachment cycles).

CONTROLS AND DISPLAY

Figure 1 - Trellix Detachment Controller Diagram

Light Indicators

- System Ready (Solid Green)
- System Ready (Blinking Green)
- Progress Indicators (Blinking Green)
- Progress Indicators (Solid Green)
- Progress Indicators (Blinking Amber)
- Low Battery (Amber)

Audio Indicators

- Single Short Beep
- Three Short Beeps
- Single Long Beep

Table 1 - Trellix Detachment Controller - Actions Table and Troubleshooting

INDICATOR	CAUSE	CONDITION / RECOMMENDED ACTIONS
No Lights	Power is Off	Delivery Pusher Not Inserted. Insert Delivery Pusher to Activate Controller.
System Ready Blinks Green	Incomplete Cable Connection (i.e. Verify Return/Ground Cable)	Return/Ground Cable Not Connected to Patient Check Cable Connections Check/Replace Grounding Needle
Single Short Beep System Ready is Solid Green	Ready to Initiate Detachment	Press Detach Button to Initiate Detachment Sequence (Must depress button for > 0.25 sec)
Progress Indicators Blink and Light Up in Sequence Over Time	Detachment in Progress	Wait
Three Short Beeps All Four (4) Progress Indicator Lights are Solid Green	Successful Detachment Detected	Verify Successful Detachment of Implant Under Fluoroscopy, then Remove Delivery Pusher.
Single Long Beep All Four (4) Progress Indicator Lights Blink Amber	Unsuccessful Detachment (Time-out)	If Ready Light is Solid Green, then Press Detach Button Again to Re-initiate Detachment
Single Long Beep All Lights Go Off	Not Ready to Initiate Detachment	Delivery Pusher Not Fully Inserted
Low Battery Light is Amber	Low Battery or Other Error	Change to New Controller

DIRECTIONS FOR USE

TRELLIX EMBOLIC COIL SYSTEM

Please refer to the Trellix Embolic Coil System Instructions for Use.

DIRECTIONS FOR USE

TRELLIX DETACHMENT CONTROLLER PREPARATION

1. Remove the Trellix Detachment Controller from its protective packaging. Place the detachment controller in the sterile field.
CAUTION: DO NOT use any power source other than the SHAPE MEMORY MEDICAL Trellix Detachment Controller to detach the SHAPE MEMORY MEDICAL Trellix Embolic Coil System. The Trellix Detachment Controller is intended to be single patient use only. DO NOT attempt to re-sterilize or otherwise re-use.
2. Insert sterile, uncoated hypodermic needle size 20 or 22 gauge (0.9mm or 0.7mm) at the patient's groin.
3. Clip the return electrode (grounding) cable (Blue cable) on the Trellix Detachment Controller onto the hypodermic needle.

DETACHMENT OF COIL

4. Once the coil has been placed into a satisfactory location, verify the hypodermic needle is in place at the patient's groin and the return electrode (grounding) blue cable connector of the Trellix Detachment Controller is connected to the hypodermic needle.
5. Verify that the RHV is firmly locked around the Trellix delivery pusher prior to insertion.
6. Prior to insertion, wipe clean and dry the exposed proximal section of the delivery pusher.
7. To insert, slide the Trellix Detachment Controller over the proximal connector of the Trellix delivery pusher firmly into the funnel section of the Trellix Detachment Controller (See Figure 2).

Figure 2 - Detachment Controller Insertion to Delivery Pusher



Figure 3 - Radiopaque Marker Alignment with Microcatheter

8. Tighten the RHV to prevent movement of the delivery pusher.
9. Verify that the distal shaft of the microcatheter is not under stress before coil detachment. Axial compression or tension forces could be stored in the microcatheter causing the tip to move during coil delivery.
CAUTION: Microcatheter tip movement could cause the aneurysm or vessel to rupture.
10. When the Trellix Detachment Controller is properly connected to the Trellix Coil System, a single short audible beep will sound and a solid green light ("System Ready") will be present on the Trellix Detachment Controller indicating the system is ready for detachment. If "System Ready" light is flashing green, check to make sure that the connection has been made and grounding cable is connected correctly.
If "System Ready" light is blinking and no audible beep was present, check and verify return electrode (grounding) cable is properly connected.
: If "Low Battery" indicator is illuminated amber, replace the Trellix Detachment Controller.
11. Prior to initiating detachment:
 - Verify visually that the continuous flush solution is infusing normally.
 - Verify under fluoroscopy, that the proximal marker of the delivery pusher is exactly distal to the proximal marker of the microcatheter with 2 tip markers (See Figure 3).
 - Verify the coil position is stable at the lesion site.
12. Press and hold momentarily the "Detach Button" on the Trellix Detachment Controller to initiate detachment. A single short audible beep will sound and the "Progress Indicators" will flash green indicating detachment. Once initiated, release the button.
13. During the detachment cycle, the "Progress Indicators" lights will progressively light up (up to a maximum duration of 120 seconds).
14. At the end of a successful detachment cycle:
 - Current flow will stop
 - Three (3) short audible beeps will sound
 - All four (4) Progress Indicator lights will be a solid green light
 - System Ready light will be off
15. If detachment is NOT successful due to a time-out (greater than 120 seconds):
 - Current flow will stop
 - One (1) long audible beep will sound
 - All four (4) Progress Indicator lights will flash amber
 - System Ready light will be a solid green

To reinitiate, repeat and press the "Detach Button" to initiate another cycle of detachment.

CAUTION: To minimize disturbance to the coil placement and delivery pusher, only press the "Detach Button" again to reinitiate. DO NOT remove or modify delivery pusher position.

16. Once the detachment of the coil has been signaled, VERIFY UNDER FLUOROSCOPY THAT THE COIL HAS DETACHED. Slowly pull back on the delivery pusher while watching fluoroscopy to make sure the coil does not move.
17. In the unlikely event the coil moves, perform the following:
- Advance the delivery pusher to re-establish the RO marker and microcatheter marker alignment.
 - Verify that the RHV is firmly locked around the Trellix delivery pusher before removing and inserting into the Trellix Detachment Controller to ensure that the coil does not move during the connection process.
 - Resume current flow by pressing and holding momentarily the "Detach Button" on the Trellix Detachment Controller.
 - Verify coil detachment under fluoroscopy as above. Repeat as necessary.
18. Once coil detachment has been detected and fluoroscopically confirmed, remove the delivery pusher from the Trellix Detachment Controller and slowly withdraw the wire from the microcatheter.
- CAUTION:** DO NOT advance the delivery pusher once the coil has been detached as it may increase risk of aneurysm or vessel perforation.
19. Repeat steps if additional coil placement is required.

TECHNICAL DESCRIPTION

Power	6V from four (4) series alkaline batteries	Unit Weight	80 g (2.8 oz)
Battery	Four (4) AAAA (1.5V) alkaline batteries, factory installed	Detachment Timeout	120 seconds
Expected Battery Life	Minimum 20 detachments	Power Source Location (IEC 60601-1)	Internally Powered Equipment
Current	1.8 mA maximum direct current (DC)	Device Classification (IEC 60601-1)	Type BF
Voltage	28 VDC maximum output	Compatibility with Flammable Anaesthetic Mixtures (IEC 60601-1)	Equipment not suitable for use in the presence of a Flammable Anaesthetic Mixture or with Oxygen or Nitrous Oxide
Operating Environment	Temperature: 10°C to 40°C Humidity: 10% to 90% RH, non-condensing		
Storage/Transportation Environment	Temperature: -18°C to 55°C Humidity: 30% to 85% RH, non-condensing	Protection against Ingress of Water (IEC 60601-1)	IPX0
Unit Size	4.2 cm (1.7 in) X 16.8 cm (6.6 in) X 2.6 cm (1.0 in)	Mode of Operation (IEC 60601-1)	Continuous Operation

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The Trellix Detachment Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the Trellix Detachment Controller should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions, CISPR 11	Group 1	The Trellix Detachment Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions, CISPR 11	Class B	
Harmonic emissions, IEC61000-3-2	Not applicable	The emissions characteristics of this equipment make it suitable for use either in industrial areas and hospitals (CISPR 11 Class A) or in a residential environment (for which CISPR 11 Class B is normally required).
Voltage fluctuations/flicker emissions, IEC 61000-3-3	Not applicable	

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Trellix Detachment Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the Trellix Detachment Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD), IEC 61000-4-2	15 kV Air Discharge (max.) 8 kV Contact Discharge	15 kV Air Discharge (max.) 8 kV Contact Discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst, IEC 61000-4-4	2 kV AC Mains 1 kV I/O Ports 100 kHz PRR	Not applicable	
Surges - AC Mains, IEC 61000-4-5	2 kV (max.) 0, 90, 180 & 270 degree phase angles	Not applicable	Not applicable
Surges - 12 VDC Power, ISO 7637-2	600 V	Not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines, IEC 61000-4-11	<ul style="list-style-type: none"> • $U_i = 0\%$, 0.5 cycle (0, 45, 90, 135, 180, 225, 270 and 315°) • $U_i = 0\%$, 1 cycle • $U_i = 70\%$, 25/30 cycles (@ 0°) • U_i, 250/300 cycles 	Not applicable	
Power frequency (50/60 Hz) magnetic field, IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_i is the AC mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Trellix Detachment Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the Trellix Detachment Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz 6Vrms (ISM bands)	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the Trellix Detachment Controller, including cables, than the recommended separation distance calculation from the equation applicable to the frequency of the transmitter. Recommended separation distance Not Applicable $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz 80% @ 1kHz AM Modulation 9 V/m to 28 V/m	3 V/m 80 MHz to 2.7 GHz 80% @ 1kHz AM Modulation 9 V/m to 28 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Trellix Detachment Controller is used exceeds the applicable RF compliance level above, the Trellix Detachment Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Trellix Detachment Controller.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE TRELIX DETACHMENT CONTROLLER

The Trellix Detachment Controller is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Trellix Detachment Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Trellix Detachment Controller as recommended below, according to the maximum output power of the communications equipment.

	Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
		150 kHz to 80 MHz $d=(3.5V1)\sqrt{P}$	80 MHz to 800 MHz $d=(3.5E1)\sqrt{P}$
0.01	not applicable	0.12	0.23
0.1	not applicable	0.38	0.73
1	not applicable	1.2	2.3
10	not applicable	3.8	7.3
100	not applicable	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT DISPOSAL INFORMATION

After the device usage, the Trellix Detachment Controller should be disposed of according to normal hospital practices (e.g., biohazardous materials should be properly handled and disposed of in appropriate containers). Infected medical devices are not regulated under the WEEE Directive but rather must be disposed of as medical waste.

WARRANTY

Shape Memory Medical Inc. warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Shape Memory Medical Inc.'s obligation under this warranty is limited to replacing or repairing at its option, this product if returned within the warranty period to Shape Memory Medical Inc. and after confirmed to be defective by the manufacturer.

Except as expressly provided in this warranty, Shape Memory Medical Inc. disclaims any representation or warranty of any kind, express or implied, including any warranty as to merchantability or fitness for a particular purpose.

See the Terms and Conditions of Sale for further information.

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